

Tranquility™

Instructions to perform the test

Always store the *Softbox* gel pack, the blood tubes and the box at room temperature. **Do not refrigerate or freeze.**

Book the courier **1 day before** the blood draw following the *Instructions to ship the sample*.

Login to www.genoma.com and **activate your kit** before shipping the sample.

Check the kit

1. Ensure that all documents and materials are present:



Test requisition form and informed consent (TRF)



Instructions to ship the sample



3 Proforma Invoice



DHL Shipment Waybill




DHL Express Flyer shipping bag



EPS box



Card to activate your kit



4 or 5 Free labels (AGxxxxxxxxxxx)



Biohazard bag and absorbent sheet



1 or 2 Blood tubes



Softbox gel pack Do not refrigerate

2. Check the expiry date on the blood tube(s)
If expired, please contact us: +41 58 1000 100 or genoma@genoma.com

Fill out the documents


3. Fill out all the fields of the Test requisition form and informed consent (TRF) and:

- write clearly on the top of the TRF the 14 digits from the tube(s) barcode (AGxxxxxxxxxxx)
- ensure your prescriber has signed the TRF
- sign on pages 1 and 2 of the TRF

4. Fill out all the fields of the Card to activate your kit and write clearly the 14 digits from the tube(s) barcode (AGxxxxxxxxxxx)

Blood draw

- Ask the qualified person for a blood draw to fill the provided tube(s) with 8 to 10 mL of venous blood per tube
- Mix by gentle inversion 8-10 times



- Stick the provided labels (AGxxxxxxxxxxx):

- one 'TUBE' label on each tube
- one 'TRF' label on the each side of the TRF
- one 'CARD' label on the *Card to activate your kit*

1

2

Front Back

3

Stick the 'TUBE' label under the word *Streck* without covering the existing barcode (GAxxxxxxxxxxx)

Now that the *Card to activate your kit* is completed, keep it for yourself

Pack and ship

- Place the labelled tube(s) into the biohazard bag which contains the absorbent sheet, and seal it carefully (Do not refrigerate)
- Place the Softbox gel pack and the biohazard bag containing the tube(s) into the EPS box and close it. Then, place it inside Tranquility box, add the labelled TRF on top of the EPS box and close the Tranquility box
- Ship the Tranquility box following the Instructions to ship the sample



The report will be issued within 4 working days on average from the time of reception in the laboratory (unless further analysis are required)

For any question you may have, please contact us: +41 58 1000 100 or genoma@genoma.com